

# PENN ENDO CREATION INSTRUMENTS

according to Dr. Syngcuk Kim



**8307**

## **Penn Endo Exam / Diagnostic Set**

9218

Mouth mirror handle PEEK, black, SIMPLE STEM

5597

Rhodium 4 - DOUBLE MIRROR (6 pcs.)

8450

Explorer DG16 / Perio Probe CP 12, PEEK black

8451

Endo Explorer, KX 1, PEEK black

4140

Pliers SEMKEN-TAYLOR straight, 15 cm

7231

Scalpel blade holder TRINOVO, blade in 7 positions adjustable

7600-BL

Steritray MINI with transverse bar



Dr. Kim

“ The perio probes and explorers are important not only for checking perio probing and related measurements but additionally to check the height of tooth. By placing the explorer on the buccal surface of the tooth, one can estimate the height of the tooth where the osteotomy must be made ”

Kohdent Roland Kohler Medizintechnik GmbH & Co. KG

Bodenseeallee 14-16 | 78333 Stockach | Germany








Tel.: +49 7771 64999-0 | Fax: +49 7771 64999-50 | info@kohler-medizintechnik.de

www.kohler-medizintechnik.de

# PENN ENDO CREATION INSTRUMENTS

according to Dr. Syngcuk Kim

## 8307 Penn Endo Exam / Diagnostic Set

	9218	Mouth mirror handle PEEK, black, SIMPLE STEM
	5597	Rhodium 4 - DOUBLE MIRROR (6 pcs.)
	8450	Explorer DG16 / Perio Probe CP 12, PEEK black
	8451	Endo Explorer, KX 1, PEEK black
	4140	Pliers SEMKEN-TAYLOR straight, 15 cm
	7231	Scalpel blade holder TRINOVO, blade in 7 positions adjustable
	7600-BL	Steritray MINI with transverse bar