

PENN ENDO CREATION INSTRUMENTS

according to Dr. Syngcuk Kim



8309 Penn Endo Plugger Set

8420	MTA Block acc. to Dr. Lee
8460	Carver / Filling instrument, PEEK black
8461	Perforation repair instrument, PEEK black
8462	Micro-Plugger, length 3 mm, 0.6 mm dia. PEEK black
8463	Micro-Plugger, length 2.5 mm, 1.0 mm dia. PEEK black
7600-BL	Steritray MINI with transverse bar



Dr. Kim

“ Plugging instruments: After placement of bioceramic materials or other filling materials, the root prepared canals must be gently condensed coronally. The Plastic Instruments above place the filling material into the root prepared canal. There are two sizes of pluggers - thin and thick – which are used to plug the materials deeply into the canal. The size of the micro-condenser tips used depends upon the size of root end prep. The plastic instrument is designed to press the root end filling material into the root end prep space. The MTA Block is a special system which allows the dentist to pick small piece of MTA or bioceramic material from the plastic block with the plastic instrument and carry to the root end prep canal! ”

PENN ENDO CREATION INSTRUMENTS

according to Dr. Syngcuk Kim

8309 Penn Endo Plugger Set

	8420	MTA Block acc. to Dr. Lee
	8460	Carver / Filling instrument, PEEK black
	8461	Perforation repair instrument, PEEK black
	8462	Micro-Plugger, length 3 mm, 0.6 mm dia. PEEK black
	8463	Micro-Plugger, length 2.5 mm, 1.0 mm dia. PEEK black
	7600-BL	Steritray MINI with transverse bar