

Granulation Tissue Instruments

acc. to Dr. Mike Sabeti, DDS
University of California, San Francisco



3305

2 mm diameter,
short and long
PEEK handle, yellow



3306

4 mm diameter,
short and long
PEEK handle, green

This surgical device is capable of efficiently and effectively removing the granulation tissue encountered in periodontal surgery, endodontic surgery and extraction socket.

The current surgical devices currently utilized in the surgical removal of granulation tissue are extremely time consuming. The granulation tissue devices have been designed to effectively remove granulation tissue. The proposed surgical device consists of a loop in the middle with sharp edges that grab and cut the granulation tissue fast and more efficient.

